

Antibodies Product Datasheet

Product Name:	CHST1 Antibody
Catalog Number:	C14920
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CHST1.
Concentration:	1 mg/ml
Gene Name:	CHST1
Protein Name:	Carbohydrate sulfotransferase 1
Other Name:	CHST1; Carbohydrate sulfotransferase 1; Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 1; GST-1; Keratan sulfate Gal-6 sulfotransferase; KS6ST; KSGal6ST; KSST
Human Swiss Prot No.:	O43916
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHST1 Polyclonal Antibody detects endogenous levels of CHST1 protein.
Storage and Stability:	-20°C/1 year