

Antibodies Product Datasheet

Product Name:	CHST6 Antibody
Catalog Number:	C14929
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHST6. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	CHST6
Protein Name:	Carbohydrate sulfotransferase 6
Other Name:	CHST6; Carbohydrate sulfotransferase 6; Corneal N-acetylglucosamine-6-O-sulfotransferase; C-GlcNAc6ST; hCGn6ST; Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 4-beta; GST4-beta; N-acetylglucosamine 6-O-sulfotransfera
Human Swiss Prot No.:	Q9GZX3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHST6 Polyclonal Antibody detects endogenous levels of CHST6 protein.
Storage and Stability:	-20°C/1 year