

Antibodies Product Datasheet

Product Name:	GalNAc4ST-1 Antibody
Catalog Number:	C14931
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHST8. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	CHST8
Protein Name:	Carbohydrate sulfotransferase 8
Other Name:	CHST8; Carbohydrate sulfotransferase 8; GalNAc-4-O-sulfotransferase 1; GalNAc-4-ST1; GalNAc4ST-1; N-acetylgalactosamine-4-O-sulfotransferase 1
Human Swiss Prot No.:	Q9H2A9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GalNAc4ST-1 Polyclonal Antibody detects endogenous levels of GalNAc4ST-1 protein.
Storage and Stability:	-20°C/1 year