

Antibodies Product Datasheet

Product Name:	CHST9 Antibody
Catalog Number:	C14932
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHST9. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	CHST9
Protein Name:	Carbohydrate sulfotransferase 9
Other Name:	CHST9; Carbohydrate sulfotransferase 9; GalNAc-4-O-sulfotransferase 2; GalNAc-4-ST2; GalNAc4ST-2; N-acetylgalactosamine-4-O-sulfotransferase 2
Human Swiss Prot No.:	Q7L1S5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHST9 Polyclonal Antibody detects endogenous levels of CHST9 protein.
Storage and Stability:	-20°C/1 year