

Antibodies Product Datasheet

Product Name:	CA VI Antibody
Catalog Number:	C14943
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CA6. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	CA6
Protein Name:	Carbonic anhydrase 6
Other Name:	CA6; Carbonic anhydrase 6; Carbonate dehydratase VI; Carbonic anhydrase VI; CA-VI; Salivary carbonic anhydrase; Secreted carbonic anhydrase
Human Swiss Prot No.:	P23280
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CA VI Polyclonal Antibody detects endogenous levels of CA VI protein.
Storage and Stability:	-20°C/1 year