

Antibodies Product Datasheet

Product Name:	CA VIII Antibody
Catalog Number:	C14946
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CA VIII.
Concentration:	1 mg/ml
Gene Name:	CA8
Protein Name:	Carbonic anhydrase-related protein
Other Name:	CA8; CALS; Carbonic anhydrase-related protein; CARP; Carbonic anhydrase VIII; CA-VIII
Human Swiss Prot No.:	P35219
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CA VIII Polyclonal Antibody detects endogenous levels of CA VIII protein.
Storage and Stability:	-20°C/1 year