

Antibodies Product Datasheet

Product Name:	CES2 Antibody
Catalog Number:	C14951
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CES2.
Concentration:	1 mg/ml
Gene Name:	CES2
Protein Name:	Cocaine esterase
Other Name:	CES2; ICE; Cocaine esterase; Carboxylesterase 2; CE-2; hCE-2; Methylumbelliferyl-acetate deacetylase 2
Human Swiss Prot No.:	O00748
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CES2 Polyclonal Antibody detects endogenous levels of CES2 protein.
Storage and Stability:	-20°C/1 year