

Antibodies Product Datasheet

Product Name:	CPN cat Antibody
Catalog Number:	C14963
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CPN1. AA range:409-458
Concentration:	1 mg/ml
Gene Name:	CPN1
Protein Name:	Carboxypeptidase N catalytic chain
Other Name:	CPN1; ACBP; Carboxypeptidase N catalytic chain; CPN; Anaphylatoxin inactivator; Arginine carboxypeptidase; Carboxypeptidase N polypeptide 1; Carboxypeptidase N small subunit; Kininase-1; Lysine carboxypeptidase; Plasma carboxypeptidase B; S
Human Swiss Prot No.:	P15169
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CPN cat Polyclonal Antibody detects endogenous levels of CPN cat protein.
Storage and Stability:	-20°C/1 year