

Antibodies Product Datasheet

Product Name:	CPN reg Antibody
Catalog Number:	C14964
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human CPN reg.
Concentration:	1 mg/ml
Gene Name:	CPN2
Protein Name:	Carboxypeptidase N subunit 2
Other Name:	CPN2; ACBP; Carboxypeptidase N subunit 2; Carboxypeptidase N 83 kDa chain; Carboxypeptidase N large subunit; Carboxypeptidase N polypeptide 2; Carboxypeptidase N regulatory subunit
Human Swiss Prot No.:	P22792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CPN reg Polyclonal Antibody detects endogenous levels of CPN reg protein.
Storage and Stability:	-20°C/1 year