

Antibodies Product Datasheet

Product Name:	Cdc42EP3 Antibody
Catalog Number:	C15024
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human Cdc42EP3.
Concentration:	1 mg/ml
Gene Name:	CDC42EP3
Protein Name:	Cdc42 effector protein 3
Other Name:	CDC42EP3; BORG2; CEP3; Cdc42 effector protein 3; Binder of Rho GTPases 2; MSE55-related Cdc42-binding protein
Human Swiss Prot No.:	Q9UKI2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc42EP3 Polyclonal Antibody detects endogenous levels of Cdc42EP3 protein.
Storage and Stability:	-20°C/1 year