

Antibodies Product Datasheet

Product Name:	CEP290 Antibody
Catalog Number:	C15036
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CEP290. AA range:771-820
Concentration:	1 mg/ml
Gene Name:	CEP290
Protein Name:	Centrosomal protein of 290 kDa
Other Name:	CEP290; BBS14; KIAA0373; NPHP6; Centrosomal protein of 290 kDa; Cep290; Bardet-Biedl syndrome 14 protein; Cancer/testis antigen 87; CT87; Nephrocystin-6; Tumor antigen se2-2
Human Swiss Prot No.:	O15078
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CEP290 Polyclonal Antibody detects endogenous levels of CEP290 protein.
Storage and Stability:	-20°C/1 year