

Antibodies Product Datasheet

Product Name:	CEP170 Antibody
Catalog Number:	C15040
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CEP170. AA range:701-750
Concentration:	1 mg/ml
Gene Name:	CEP170
Protein Name:	Centrosomal protein of 170 kDa
Other Name:	CEP170; FAM68A; KAB; KIAA0470; Centrosomal protein of 170 kDa; Cep170; KARP-1-binding protein; KARP1-binding protein
Human Swiss Prot No.:	Q5SW79
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CEP170 Polyclonal Antibody detects endogenous levels of CEP170 protein.
Storage and Stability:	-20°C/1 year