

Antibodies Product Datasheet

Product Name:	CEP57 Antibody
Catalog Number:	C15044
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CEP57. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	CEP57
Protein Name:	Centrosomal protein of 57 kDa
Other Name:	CEP57; KIAA0092; TSP57; Centrosomal protein of 57 kDa; Cep57; FGF2-interacting protein; Testis-specific protein 57; Translokin
Human Swiss Prot No.:	Q86XR8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CEP57 Polyclonal Antibody detects endogenous levels of CEP57 protein.
Storage and Stability:	-20°C/1 year