

Antibodies Product Datasheet

Product Name:	CLCC1 Antibody
Catalog Number:	C15076
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLCC1. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	CLCC1
Protein Name:	Chloride channel CLIC-like protein 1
Other Name:	CLCC1; KIAA0761; MCLC; Chloride channel CLIC-like protein 1; Mid-1-related chloride channel protein 1
Human Swiss Prot No.:	Q96S66
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLCC1 Polyclonal Antibody detects endogenous levels of CLCC1 protein.
Storage and Stability:	-20°C/1 year