

Antibodies Product Datasheet

Product Name:	CLC-4 Antibody
Catalog Number:	C15079
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLCN4. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	CLCN4
Protein Name:	H(+)/Cl(-) exchange transporter 4
Other Name:	CLCN4; H(+)/Cl(-) exchange transporter 4; Chloride channel protein 4; CLC-4; Chloride transporter CLC-4
Human Swiss Prot No.:	P51793
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLC-4 Polyclonal Antibody detects endogenous levels of CLC-4 protein.
Storage and Stability:	-20°C/1 year