

Antibodies Product Datasheet

Product Name:	CLC-6 Antibody
Catalog Number:	C15081
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLCN6. AA range:611-660
Concentration:	1 mg/ml
Gene Name:	CLCN6
Protein Name:	Chloride transport protein 6
Other Name:	CLCN6; KIAA0046; Chloride transport protein 6; Chloride channel protein 6; CIC-6
Human Swiss Prot No.:	P51797
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLC-6 Polyclonal Antibody detects endogenous levels of CLC-6 protein.
Storage and Stability:	-20°C/1 year