

Antibodies Product Datasheet

Product Name:	CLC-7 Antibody
Catalog Number:	C15082
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLCN7. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	CLCN7
Protein Name:	H(+)/Cl(-) exchange transporter 7
Other Name:	CLCN7; H(+)/Cl(-) exchange transporter 7; Chloride channel 7 alpha subunit; Chloride channel protein 7; CLC-7
Human Swiss Prot No.:	P51798
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLC-7 Polyclonal Antibody detects endogenous levels of CLC-7 protein.
Storage and Stability:	-20°C/1 year