

Antibodies Product Datasheet

Product Name:	CLC-KA Antibody
Catalog Number:	C15083
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLCNKA. AA range:581-630
Concentration:	1 mg/ml
Gene Name:	CLCNKA
Protein Name:	Chloride channel protein CLC-Ka
Other Name:	CLCNKA; Chloride channel protein CLC-Ka; Chloride channel Ka; CLC-K1
Human Swiss Prot No.:	P51800
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLC-KA Polyclonal Antibody detects endogenous levels of CLC-KA protein.
Storage and Stability:	-20°C/1 year