

Antibodies Product Datasheet

Product Name:	CLIC6 Antibody
Catalog Number:	C15086
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLIC6. AA range:411-460
Concentration:	1 mg/ml
Gene Name:	CLIC6
Protein Name:	Chloride intracellular channel protein 6
Other Name:	CLIC6; CLIC1L; Chloride intracellular channel protein 6; Parchorin
Human Swiss Prot No.:	Q96NY7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLIC6 Polyclonal Antibody detects endogenous levels of CLIC6 protein.
Storage and Stability:	-20°C/1 year