

Antibodies Product Datasheet

Product Name:	CSGlcA-T Antibody
Catalog Number:	C15095
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSGLCAT. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CHPF2
Protein Name:	Chondroitin sulfate glucuronyltransferase
Other Name:	CHPF2; CHSY3; CSGLCAT; KIAA1402; Chondroitin sulfate glucuronyltransferase; CSGlcA-T; Chondroitin glucuronyltransferase; Chondroitin polymerizing factor 2; ChPF-2; Chondroitin synthase 3; ChSy-3; N-acetylgalactosaminyl-proteoglycan 3-beta-g
Human Swiss Prot No.:	Q9P2E5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CSGlcA-T Polyclonal Antibody detects endogenous levels of CSGlcA-T protein.
Storage and Stability:	-20°C/1 year