

Antibodies Product Datasheet

Product Name:	CHSY1 Antibody
Catalog Number:	C15098
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHSY1. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	CHSY1
Protein Name:	Chondroitin sulfate synthase 1
Other Name:	CHSY1; CHSY; CSS1; KIAA0990; Chondroitin sulfate synthase 1; Chondroitin glucuronyltransferase 1; Chondroitin synthase 1; ChSy-1; Glucuronosyl-N-acetylgalactosaminyl-proteoglycan 4-beta-N-acetylgalactosaminyltransferase 1; N-acetylgalactosa
Human Swiss Prot No.:	Q86X52
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHSY1 Polyclonal Antibody detects endogenous levels of CHSY1 protein.
Storage and Stability:	-20°C/1 year