

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                  |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CHSY2 Antibody                                                                                                                                                                                                                                   |
| <b>Catalog Number:</b>        | C15099                                                                                                                                                                                                                                           |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                                                                                                |
| <b>Applications:</b>          | WB, IHC-p, IF/ICC, ELISA                                                                                                                                                                                                                         |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.                                                                                     |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                                                                                 |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                            |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human CHSY2. AA range:631-680                                                                                                                                                |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                          |
| <b>Gene Name:</b>             | CHPF                                                                                                                                                                                                                                             |
| <b>Protein Name:</b>          | Chondroitin sulfate synthase 2                                                                                                                                                                                                                   |
| <b>Other Name:</b>            | CHPF; CSS2; Chondroitin sulfate synthase 2; Chondroitin glucuronyltransferase 2; Chondroitin-polymerizing factor; ChPF; Glucuronosyl-N-acetylgalactosaminyl-proteoglycan 4-beta-N-acetylgalactosaminyltransferase II; N-acetylgalactosaminyl-pro |
| <b>Human Swiss Prot No.:</b>  | Q8IZ52                                                                                                                                                                                                                                           |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                          |
| <b>Specificity:</b>           | CHSY2 Polyclonal Antibody detects endogenous levels of CHSY2 protein.                                                                                                                                                                            |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                     |