

Antibodies Product Datasheet

Product Name:	CBX2 Antibody
Catalog Number:	C15106
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human CBX2.
Concentration:	1 mg/ml
Gene Name:	CBX2
Protein Name:	Chromobox protein homolog 2
Other Name:	CBX2; Chromobox protein homolog 2
Human Swiss Prot No.:	Q14781
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CBX2 Polyclonal Antibody detects endogenous levels of CBX2 protein.
Storage and Stability:	-20°C/1 year