

Antibodies Product Datasheet

Product Name:	CDYL2 Antibody
Catalog Number:	C15112
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunocytochemistry: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDYL2. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	CDYL2
Protein Name:	Chromodomain Y-like protein 2
Other Name:	CDYL2; Chromodomain Y-like protein 2; CDY-like 2
Human Swiss Prot No.:	Q8N8U2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CDYL2 Polyclonal Antibody detects endogenous levels of CDYL2 protein.
Storage and Stability:	-20°C/1 year