

Antibodies Product Datasheet

Product Name:	CMTM1 Antibody
Catalog Number:	C15131
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CKLF1. AA range:120-169
Concentration:	1 mg/ml
Gene Name:	CMTM1
Protein Name:	CKLF-like MARVEL transmembrane domain-containing protein 1
Other Name:	CMTM1; CKLF1; CKLF-like MARVEL transmembrane domain-containing protein 1; Chemokine-like factor superfamily member 1
Human Swiss Prot No.:	Q8IZ96
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CMTM1 Polyclonal Antibody detects endogenous levels of CMTM1 protein.
Storage and Stability:	-20°C/1 year