

Antibodies Product Datasheet

Product Name:	CMTM3 Antibody
Catalog Number:	C15133
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CKLF3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CMTM3
Protein Name:	CKLF-like MARVEL transmembrane domain-containing protein 3
Other Name:	CMTM3; CKLF3; CKLF-like MARVEL transmembrane domain-containing protein 3; Chemokine-like factor superfamily member 3
Human Swiss Prot No.:	Q96MX0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CMTM3 Polyclonal Antibody detects endogenous levels of CMTM3 protein.
Storage and Stability:	-20°C/1 year