

Antibodies Product Datasheet

Product Name:	Claudin-8 Antibody
Catalog Number:	C15163
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLDN8. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	CLDN8
Protein Name:	Claudin-8
Other Name:	CLDN8; Claudin-8
Human Swiss Prot No.:	P56748
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Claudin-8 Polyclonal Antibody detects endogenous levels of Claudin-8 protein.
Storage and Stability:	-20°C/1 year