

Antibodies Product Datasheet

Product Name:	CstF-64T Antibody
Catalog Number:	C15167
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSTF2T. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	CSTF2T
Protein Name:	Cleavage stimulation factor subunit 2 tau variant
Other Name:	CSTF2T; KIAA0689; Cleavage stimulation factor subunit 2 tau variant; CF-1 64 kDa subunit tau variant; Cleavage stimulation factor 64 kDa subunit tau variant; CSTF 64 kDa subunit tau variant; TauCstF-64
Human Swiss Prot No.:	Q9H0L4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CstF-64T Polyclonal Antibody detects endogenous levels of CstF-64T protein.
Storage and Stability:	-20°C/1 year