

Antibodies Product Datasheet

Product Name:	CD75 Antibody
Catalog Number:	C15170
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ST6GAL1. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	ST6GAL1
Protein Name:	Beta-galactoside alpha-2,6-sialyltransferase 1
Other Name:	ST6GAL1; SIAT1; Beta-galactoside alpha-2, 6-sialyltransferase 1; Alpha 2,6-ST 1; B-cell antigen CD75; CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,6-sialyltransferase 1; ST6Gal I; ST6GalII; Sialyltransferase 1
Human Swiss Prot No.:	P15907
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD75 Polyclonal Antibody detects endogenous levels of CD75 protein.
Storage and Stability:	-20°C/1 year