

Antibodies Product Datasheet

Product Name:	RCL Antibody
Catalog Number:	C15173
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RCL. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	RCL
Protein Name:	Deoxyribonucleoside 5'-monophosphate N-glycosidase
Other Name:	RCL; C6orf108; Deoxyribonucleoside 5'-monophosphate N-glycosidase; c-Myc-responsive protein Rcl
Human Swiss Prot No.:	O43598
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RCL Polyclonal Antibody detects endogenous levels of RCL protein.
Storage and Stability:	-20°C/1 year