

Antibodies Product Datasheet

Product Name:	Factor XIII B Antibody
Catalog Number:	C15182
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human F13B. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	F13B
Protein Name:	Coagulation factor XIII B chain
Other Name:	F13B; Coagulation factor XIII B chain; Fibrin-stabilizing factor B subunit; Protein-glutamine gamma-glutamyltransferase B chain; Transglutaminase B chain
Human Swiss Prot No.:	P05160
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Factor XIII B Polyclonal Antibody detects endogenous levels of Factor XIII B protein.
Storage and Stability:	-20°C/1 year