

Antibodies Product Datasheet

Product Name:	CCDC109A Antibody
Catalog Number:	C15187
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CCDC109A.
Concentration:	1 mg/ml
Gene Name:	MCU
Protein Name:	Calcium uniporter protein mitochondrial
Other Name:	MCU; C10orf42; CCDC109A; Calcium uniporter protein; mitochondrial; Coiled-coil domain-containing protein 109A
Human Swiss Prot No.:	Q8NE86
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCDC109A Polyclonal Antibody detects endogenous levels of CCDC109A protein.
Storage and Stability:	-20°C/1 year