

Antibodies Product Datasheet

Product Name:	CCDC4 Antibody
Catalog Number:	C15190
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BEND4. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	BEND4
Protein Name:	BEN domain-containing protein 4
Other Name:	BEND4; CCDC4; BEN domain-containing protein 4; Coiled-coil domain-containing protein 4
Human Swiss Prot No.:	Q6ZU67
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCDC4 Polyclonal Antibody detects endogenous levels of CCDC4 protein.
Storage and Stability:	-20°C/1 year