

Antibodies Product Datasheet

Product Name:	CCDC45 Antibody
Catalog Number:	C15191
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCDC45. AA range:591-640
Concentration:	1 mg/ml
Gene Name:	CEP95
Protein Name:	Centrosomal protein of 95 kDa
Other Name:	CEP95; CCDC45; CEP45; Centrosomal protein of 95 kDa; Coiled-coil domain-containing protein 45
Human Swiss Prot No.:	Q96GE4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCDC45 Polyclonal Antibody detects endogenous levels of CCDC45 protein.
Storage and Stability:	-20°C/1 year