

Antibodies Product Datasheet

Product Name:	CCDC99 Antibody
Catalog Number:	C15193
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCDC99. AA range:541-590
Concentration:	1 mg/ml
Gene Name:	CCDC99
Protein Name:	Protein Spindly
Other Name:	CCDC99; Protein Spindly; hSpindly; Arsenite-related gene 1 protein; Coiled-coil domain-containing protein 99; Rhabdomyosarcoma antigen MU-RMS-40.4A
Human Swiss Prot No.:	Q96EA4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCDC99 Polyclonal Antibody detects endogenous levels of CCDC99 protein.
Storage and Stability:	-20°C/1 year