

Antibodies Product Datasheet

Product Name:	C1QBP Antibody
Catalog Number:	C15217
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C1QBP. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	C1QBP
Protein Name:	Complement component 1 Q subcomponent-binding protein mitochondrial
Other Name:	C1QBP; GC1QBP; HABP1; SF2P32; Complement component 1 Q subcomponent-binding protein; mitochondrial; GC1q-R protein; Glycoprotein gC1qBP; C1qBP; Hyaluronan-binding protein 1; Mitochondrial matrix protein p32; p33
Human Swiss Prot No.:	Q07021
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C1QBP Polyclonal Antibody detects endogenous levels of C1QBP protein.
Storage and Stability:	-20°C/1 year