

Antibodies Product Datasheet

Product Name:	hCAP-H Antibody
Catalog Number:	C15233
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NCAPH. AA range:441-490
Concentration:	1 mg/ml
Gene Name:	NCAPH
Protein Name:	Condensin complex subunit 2
Other Name:	NCAPH; BRRN; BRRN1; CAPH; KIAA0074; Condensin complex subunit 2; Barren homolog protein 1; Chromosome-associated protein H; hCAP-H; Non-SMC condensin I complex subunit H; XCAP-H homolog
Human Swiss Prot No.:	Q15003
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	hCAP-H Polyclonal Antibody detects endogenous levels of hCAP-H protein.
Storage and Stability:	-20°C/1 year