

Antibodies Product Datasheet

Product Name:	Contactin 2 Antibody
Catalog Number:	C15238
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNTN2. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	CNTN2
Protein Name:	Contactin-2
Other Name:	CNTN2; AXT; TAG1; TAX1; Contactin-2; Axonal glycoprotein TAG-1; Axonin-1; Transient axonal glycoprotein 1; TAX-1
Human Swiss Prot No.:	Q02246
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Contactin 2 Polyclonal Antibody detects endogenous levels of Contactin 2 protein.
Storage and Stability:	-20°C/1 year