

Antibodies Product Datasheet

Product Name:	CASPR3 Antibody
Catalog Number:	C15245
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human CASPR3.
Concentration:	1 mg/ml
Gene Name:	CNTNAP3
Protein Name:	Contactin-associated protein-like 3
Other Name:	CNTNAP3; CASPR3; KIAA1714; Contactin-associated protein-like 3; Cell recognition molecule Caspr3
Human Swiss Prot No.:	Q9BZ76
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CASPR3 Polyclonal Antibody detects endogenous levels of CASPR3 protein.
Storage and Stability:	-20°C/1 year