

Antibodies Product Datasheet

Product Name:	CAR Antibody
Catalog Number:	C15257
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CXADR. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CXADR
Protein Name:	Coxsackievirus and adenovirus receptor
Other Name:	CXADR; CAR; Coxsackievirus and adenovirus receptor; CAR; hCAR; CVB3-binding protein; Coxsackievirus B-adenovirus receptor; HCVADR
Human Swiss Prot No.:	P78310
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CAR Polyclonal Antibody detects endogenous levels of CAR protein.
Storage and Stability:	-20°C/1 year