

Antibodies Product Datasheet

Product Name:	Cubilin Antibody
Catalog Number:	C15264
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human Cubilin.
Concentration:	1 mg/ml
Gene Name:	CUBN
Protein Name:	Cubilin
Other Name:	CUBN; IFCR; Cubilin; 460 kDa receptor; Intestinal intrinsic factor receptor; Intrinsic factor-cobalamin receptor; Intrinsic factor-vitamin B12 receptor
Human Swiss Prot No.:	O60494
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cubilin Polyclonal Antibody detects endogenous levels of Cubilin protein.
Storage and Stability:	-20°C/1 year