

Antibodies Product Datasheet

Product Name:	Cystatin 11 Antibody
Catalog Number:	C15270
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CST11. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	CST11
Protein Name:	Cystatin-11
Other Name:	CST11; CST8L; Cystatin-11
Human Swiss Prot No.:	Q9H112
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cystatin 11 Polyclonal Antibody detects endogenous levels of Cystatin 11 protein.
Storage and Stability:	-20°C/1 year