

Antibodies Product Datasheet

Product Name:	Cystatin C Antibody
Catalog Number:	C15276
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Cystatin C . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	CST3
Protein Name:	Cystatin-C
Other Name:	CST3; Cystatin-C; Cystatin-3; Gamma-trace; Neuroendocrine basic polypeptide; Post-gamma-globulin
Human Swiss Prot No.:	P01034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cystatin C Polyclonal Antibody detects endogenous levels of Cystatin C protein.
Storage and Stability:	-20°C/1 year