

Antibodies Product Datasheet

Product Name:	DOCK 180 Antibody
Catalog Number:	C15315
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DOCK1. AA range:1661-1710
Concentration:	1 mg/ml
Gene Name:	DOCK1
Protein Name:	Dedicator of cytokinesis protein 1
Other Name:	DOCK1; Dedicator of cytokinesis protein 1; 180 kDa protein downstream of CRK; DOCK180
Human Swiss Prot No.:	Q14185
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DOCK 180 Polyclonal Antibody detects endogenous levels of DOCK 180 protein.
Storage and Stability:	-20°C/1 year