

Antibodies Product Datasheet

Product Name:	DHRS4 Antibody
Catalog Number:	C15323
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DHRS4. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	DHRS4
Protein Name:	Dehydrogenase/reductase SDR family member 4
Other Name:	DHRS4; Dehydrogenase/reductase SDR family member 4; NADPH-dependent carbonyl reductase/NADP-retinol dehydrogenase; CR; PHCR; NADPH-dependent retinol dehydrogenase/reductase; NRDR; humNRDR; Peroxisomal short-chain alcohol dehydrogenase; PSCD
Human Swiss Prot No.:	Q9BTZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DHRS4 Polyclonal Antibody detects endogenous levels of DHRS4 protein.
Storage and Stability:	-20°C/1 year