

Antibodies Product Datasheet

Product Name:	17 β -HSD11 Antibody
Catalog Number:	C15325
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DHRS8. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	HSD17B11
Protein Name:	Estradiol 17-beta-dehydrogenase 11
Other Name:	HSD17B11; DHRS8; PAN1B; PSEC0029; Estradiol 17-beta-dehydrogenase 11; 17-beta-hydroxysteroid dehydrogenase 11; 17-beta-HSD 11; 17bHSD11; 17betaHSD11; 17-beta-hydroxysteroid dehydrogenase XI; 17-beta-HSD XI; 17betaHSDXI; Cutaneous T-cell lym Q8NBQ5
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	17 β -HSD11 Polyclonal Antibody detects endogenous levels of 17 β -HSD11 protein.
Storage and Stability:	-20°C/1 year