

Antibodies Product Datasheet

Product Name:	Dsg2 Antibody
Catalog Number:	C15333
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DSG2. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	DSG2
Protein Name:	Desmoglein-2
Other Name:	DSG2; CDHF5; Desmoglein-2; Cadherin family member 5; HDGC
Human Swiss Prot No.:	Q14126
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dsg2 Polyclonal Antibody detects endogenous levels of Dsg2 protein.
Storage and Stability:	-20°C/1 year