

Antibodies Product Datasheet

Product Name:	DGAT2L6 Antibody
Catalog Number:	C15350
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGAT2L6. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	DGAT2L6
Protein Name:	Diacylglycerol O-acyltransferase 2-like protein 6
Other Name:	DGAT2L6; DC3; Diacylglycerol O-acyltransferase 2-like protein 6; Diacylglycerol O-acyltransferase candidate 3; hDC3
Human Swiss Prot No.:	Q6ZPD8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGAT2L6 Polyclonal Antibody detects endogenous levels of DGAT2L6 protein.
Storage and Stability:	-20°C/1 year