

Antibodies Product Datasheet

Product Name:	CRMP-3 Antibody
Catalog Number:	C15357
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DPYSL4. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	DPYSL4
Protein Name:	Dihydropyrimidinase-related protein 4
Other Name:	DPYSL4; CRMP3; ULIP4; Dihydropyrimidinase-related protein 4; DRP-4; Collapsin response mediator protein 3; CRMP-3; UNC33-like phosphoprotein 4; ULIP-4
Human Swiss Prot No.:	O14531
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRMP-3 Polyclonal Antibody detects endogenous levels of CRMP-3 protein.
Storage and Stability:	-20°C/1 year