

Antibodies Product Datasheet

Product Name:	HURP Antibody
Catalog Number:	C15367
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DLGAP5. AA range:791-840
Concentration:	1 mg/ml
Gene Name:	DLGAP5
Protein Name:	Disks large-associated protein 5
Other Name:	DLGAP5; DLG7; KIAA0008; Disks large-associated protein 5; DAP-5; Discs large homolog 7; Disks large-associated protein DLG7; Hepatoma up-regulated protein; HURP
Human Swiss Prot No.:	Q15398
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HURP Polyclonal Antibody detects endogenous levels of HURP protein.
Storage and Stability:	-20°C/1 year